curea medical GmbH
Münsterstraße 61-65
48565 Steinfurt
info@curea-medical.com
+49 36071 9009500

STERILE EO

MD

CE 0044

curea medical

Innovation in Care

curea P2 active

Gebrauchsanweisung - Instructions for use - Manuel d'utilisation - Gebruiksaanwijzing - Instrucciones de uso - Istruzioni per l'uso - Instruções de utilização - Brugsanvisninger - Käyttöohjeet - Bruksanvisningar - Instruksjoner for bruk - Instrukcja użytkowania - Használati utasítás - Upute za korištenje - Návod k použití - Návod na použitie - Navodila za uporabo - Uputstvo za upotrebu - Instrucțiuni de utilizare - Kasutusjuhend - Lietošanas instrukcija - Naudojimo instrukcijos - Инструкции за употреба - Οδηγίες χρήσης - Инструкции по применению - Kullanım Talimatları - 使用说明 - مادختس الالامى لعت

Inhaltsverzeichnis/ Table of content

| | |
|-------------|----|
| DEUTSCH | 4 |
| ENGLISH | 6 |
| FRANÇAIS | 8 |
| NEDERLANDS | 10 |
| ESPAÑOL | 12 |
| ITALIANO | 14 |
| PORTUGUÊS | 16 |
| DANSK | 18 |
| SUOMALAINEN | 20 |
| SVENSKA | 22 |
| NORSK | 24 |
| POLSKI | 26 |
| MAGYAR | 28 |
| HRVATSKI | 30 |
| ČESKY | 32 |
| SLOVENSKÝ | 34 |
| SLOVENSKO | 36 |
| SRPSKI | 38 |
| ROMÂNĂ | 40 |
| EESTLANE | 42 |
| LATVIETIS | 44 |
| LIETUVIS | 46 |
| БЪЛГАРСКИ | 48 |
| ΕΛΛΗΝΙΚΑ | 50 |
| РУССКИЙ | 52 |
| TÜRK | 54 |
| 中文 | 56 |
| عربي | 58 |
| SYMBOLS | 60 |

DEUTSCH

1) Beschreibung

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ENGLISH

1) Description

curea P2 active is a sterile wound dressing based on polyacrylate. The wound-facing surface comprises a three-dimensional sheet made of polyethylene, the reverse side is made of breathable yet fluid-impermeable BTBS made of polyethylene and polypropylene. The core is a mixture of cotton, cellulose and sodium polyacrylate and a layer of activated carbon. Bacteria, toxins and odour molecules are securely bound in the core material. In use, the sodium polyacrylate forms a gel. This causes the exudate to be locked into the core without dripping or leaking. The cellulose ensures that a moist wound environment is maintained (hydrobalance). Bacteria, toxins and odour molecules are securely bound in the core material. Whole blood is safely bound in the wound dressing. curea P2 active also ensures reliable exudate absorption within the framework of compression therapy. The device is sterilised with ethylene oxide and contains neither pharmacologically active substances nor latex or phthalates. The wound dressing's specified performance is 1.5 g of exudate per cm² of core surface. The dimensions of the core surface are given on the reverse of the folded box. The overall performance is the product of the width x the length x 1.5.

2) Indication

Sterile coverage of wounds and absorption of wound exudate. Wounds with secondary healing and moderate to heavy exudation. Chronic, acute and iatrogenic wounds.

3) Contraindication

- Dry wounds as this can lead to the stagnation of wound healing or cause the wound dressing to adhere to the wound bed.
- Do not use on mucous membranes, eyes or exposed bones or tendons as these body parts may be permanently damaged as a result.
- For external application only: Care must be taken with wound pockets as the device may expand during use resulting in pain.
- Known over-sensitivity to the ingredients (acrylic acid).
- Heavily bleeding wounds since the device does not function as a haemostatic agent (see 5).

4) Application

The product is intended for use by medical or nursing professionals. Laypersons should use this product only after thorough instruction by appropriate qualified personnel. Ensure that the underlying disease is treated and apply the devices under medical supervision and in line with the indication. The time intervals for changing the dressing and inspecting the wound are determined by the attending doctor based on the specific wound situation. We recommend changing the dressing at the latest after seven days and before the dressing becomes excessively full. Depending on the infection status, we recommend regular, ideally daily inspection of the area surrounding the wound and the dressing. In the case of very heavy exudation, we recommend selecting a larger wound dressing than actually required to cover the wound surface. Ensure that there is good contact between the wound dressing and wound bed. Please remember that protein-rich, viscous exudate can block the surface of the wound dressing. In such cases, the dressing must be changed at shorter intervals. In unfavourable conditions, adhesive wound dressing can be rinsed with sterile sodium chloride solution or Ringer's solution and removed. Ensure that there is a sufficient amount of space to fix the wound dressing securely. Use standard elastic fixation materials to secure the dressing in place. When fixing the dressing, take care to prevent tension blisters or skin damage to already damaged skin (e.g., parchment skin)

5) Residual risk

As a rule, application of the devices if used as intended and subject to compliance with the instructions for use and the established treatment standards is very low risk, nevertheless the following residual risks do exist:

- Given that the wound dressing swells up to 1.5 cm, fixation bandage and self-adhesive wound dressing can cause tension on the skin (tension blisters). Particular care must be taken with patients with known skin damage or age-induced atrophy ("parchment skin").
- Excessively full wound dressings can cause maceration of the skin.
- Very full wound dressings can impair the circulation as a result of pressure and their weight.
- Infections can go unnoticed due to long application times. Check the area surrounding the wound on a regular basis for signs of infection such as redness, increased skin temperature and increased pain.
- You must be very cautious with wounds which are bleeding heavily: depending on their size, the wound dressings can absorb larger amounts of blood and thus cause a critical health condition.

6) Side effects/obligation to report

There are no known physiological side-effects. With combinations of liquid or paste-like therapeutic agents these can be absorbed by the wound dressing. The effectiveness of these devices and the wound dressings will be reduced accordingly. The application of the devices is not restricted to any specific population group. Should side effects occur in connection with the application of this device, please report these informally by post, telephone or e-mail to info@curea-medical.de.

If, in the unlikely event that the side-effects are so serious that one of the following consequences occurred directly or indirectly, could have occurred or could occur in the future:

- death of a patient, user or another person,
- temporary or permanent serious deterioration in the health of the patient, user or another person,
- a serious risk for public health,

please report these immediately to the competent supervisory authority and the manufacturer.

7) Warnings

- Do not tear or cut the wound dressing and this could cause sodium polyacrylate to escape and soil the patient. Should material leak out, clean the wound with sterile wound irrigation solution, there is no risk to either the patient or user.
- Do not use damaged, incorrectly stored or expired devices in order to avoid possible infection risks, dispose of affected devices as per the valid regulations.
- The wound dressings are intended for single use. Do not reuse in order to rule out the risk of possible infections.
- Do not resterilise in order to rule out the risk of possible infections or damage to the wound dressings.

8) Storage

Protect against moisture and direct sunlight and note the limits for humidity and temperature.

9) Disposal

The wound dressing can be disposed of as normal clinical waste or domestic waste. The local regulations must be noted.

10) Advantage/Disadvantage

Benefit: Odour reduction and binding of toxins and bacteria. Reduced risk of adherence to the wound bed.

Disadvantage: No cleaning effect due to adhesion of the wound coatings.

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